SERVOX® Basisplatte

DE

SERVOX® baseplate

EN

SERVOX® Piastra base

IT

SERVOX® podložna pločica

HR

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NL

Plaster SERVOX®

PL

Orbisana Healthcare GmbH

Biberweg 24-26 | 53842 Troisdorf | Germany

Tel. +49 (0)2241 9322-0 | Fax +49 (0)2241 9322-277

info@orbisana.de | www.orbisana.de

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SERVOX® Basisplatte

Verwendete Symbole 30-32

1. Vorwort

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EN

1. Foreword

This instruction manual is intended for patients, relatives and medical staff. SERVOX® baseplates are suitable for use at home and in a hospital setting. The patient can use the product independently (without the support of relatives or carers), but must be given instructions on its use by their doctor or by medical staff beforehand.

The SERVOX® baseplate must be used in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine.

Note

Read this instruction manual carefully. It contains important information on the safe and proper use of the product. Keep this instruction manual for future reference.

SERVOX® baseplates are assigned a batch number for the purpose of clear traceability. The batch number is indicated on the packaging next to "LOT". If you have any questions or problems with SERVOX® baseplates, please contact your field representative or Orbisana Healthcare GmbH, quoting the product's batch number.

1.1 Intended purpose

The SERVOX® baseplate is glued externally to the skin surrounding an artificially created breathing hole (tracheostoma) to permit the fitting of accessories with a 22 mm connector.

Note

SERVOX® baseplates are single-use products and must not be reused. If reused, there is a risk of infection.

26347/Version 2024-08-23

1.2 Indication

- Tracheotomised or larygectomised patients using aids with a 22 mm connection.

1.3 Contraindication

- Fresh bleeding wounds
- Chronic ulcerating wounds
- Allergies to the material used
- Use during radiotherapy or within the first two weeks after radiotherapy within the tracheostoma region.

1.4 Side effects

- Possible reduction of the size of the tracheostoma opening due to the size of the holes in the baseplate.
- Skin irritation or allergic reactions to the glue used. In this case, discontinue use of the SERVOX® baseplate and seek your doctor's advice.

2. Safety instructions

Safe and proper use of the SERVOX® baseplate is important for your health. Ignoring these safety instructions can have significant effects on your health or even cause life-threatening injuries.

In the following, warnings on the safe use of the product are indicated with the  symbol.

Danger: The word "Danger" indicates a likely situation that can result in grave injury or even death if not avoided.

Caution: The word "Caution" indicates a situation that can result in minor or moderate injury if not avoided.

Please observe the following safety instructions:

Danger

Risk of suffocation: Laryngectomised patients are unable to breathe if the tracheostoma is closed. They may experience respiratory distress or suffocate, injure themselves or die.

- Do not use speech valves or caps in laryngectomised patients that are intended for tracheotomised patients!

Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.

3. Description

The SERVOX® baseplate is a transparent and elastic self-adhesive plater with an integral 22 mm attachment ring over a central hole for the tracheostoma. The hypoallergenic and latex-free material of the plaster is coated with a medical glue, proving the tracheostoma with an all-round seal. The adhesive surface is protected with a protective foil that must be removed prior to use.

Each SERVOX® baseplate is individually wrapped and has to be replaced or removed and then disposed of latest after 24 hours or earlier if required.

The baseplate's 22 mm attachment ring can be used to fit HME cassettes (heat and moisture exchanger) or tracheostoma valves as well as other products with a 22 mm connector such as speech valves (only tracheotomised and not laryngectomised patients) or tracheal cannulas suitable for connection to a baseplate.

Various shapes of baseplates permit adaptation to the individual situation.

71 x 76 mm

SERVOX® baseplate *flex*
Stays in place even during physical activity

73 x 89 mm

SERVOX® baseplate *flex plus*
With a slightly enlarged adhesive surface

88 x 108 mm

SERVOX® baseplate *flex fit*
With a much enlarged adhesive surface

3.1 Combinable products

- Accessories with a 22 mm connection, such as speech valves (only tracheotomised and not laryngectomised patients), tracheostoma valves or HME cassettes
- SERVOX® HME Modul+
- SERVOX® 22/15 Adapter
- SERVOX® closure cap (only tracheotomised and not laryngectomised patients)
- SERVOX® Soft Tube 22 *connect*

4. Application

Caution

Risk of injury/infection: Any damage to the baseplate or its packaging may lead to injury or infection.

- Do not use the baseplate if there is even the slightest contamination to the baseplate or its packaging!

Risk of infection: Glue residue, foreign bodies or liquid that ingress into the tracheostoma during the removal of the baseplate and the

cleaning of the skin around the tracheostoma can be aspirated or cause infections.

- Proceed with caution when removing or attaching the baseplate or preparing the skin around the tracheostoma!

Note

- The SERVOX® baseplate has to be replaced or removed latest after 24 hours or earlier if required.
- If you have any questions or experience complications, please contact your doctor.
- Special tracheostoma care products (e.g. SERVOX® 3-phase-care) can be used for the cleaning and preparatory degreasing of the skin around the tracheostoma. Always follow the manufacturer's instructions for use.

4.1 Attaching the SERVOX® baseplate

- Wash your hands thoroughly and disinfect as appropriate prior to handling the baseplate.
- Check the baseplate and its packaging for any damage prior to use.
- Clean the skin around the tracheostoma with lukewarm water and a mild soap or a special tracheostoma cleaning product, and dry carefully. The skin around the tracheostoma should be as free from grease and oil as possible. If necessary, a specialised product can be used to degrease and prepare the skin. Always make sure that no water or other substances get into the trachea.
- Take the baseplate out of the packaging and remove the protective foil from the adhesive surface. Make sure that the baseplate does not fold over or that the adhesive surface gets contaminated.
- Slightly stretch the skin around the tracheostoma and stick the baseplate with the opening directly over the tracheostoma to the skin and press down. Avoid any air bubbles between the skin and the adhesive surface.
- Accessories can then be fitted to the 22 mm attachment ring.

4.2 Changing the SERVOX® baseplate

- To remove the SERVOX® baseplate, grip the tab on the side, hold down the surrounding skin and carefully peel the baseplate off the skin.
- Any glue residue on the skin can be carefully removed with an adhesive dissolver, making sure that neither glue residue nor liquid get into the tracheostoma.
- Prior to attaching a new SERVOX® baseplate, clean and, if required, pre-treat the skin.

5. Storage, transport and disposal

Store or transport the SERVOX® baseplates in a clean, dry environment that is free of dust and protected from light at temperatures between 10°C and +25°C. Protect from moisture, direct sunlight and heat. Dispose after use in domestic waste. The packaging can be recycled.

6. Technical data / Order information

Order information

Description	Size (mm)	PU	Order no.
SERVOX® baseplate <i>flex</i> *	71 x 76	30	25804
SERVOX® baseplate <i>flex</i> *	71 x 76	10	26104
SERVOX® baseplate <i>flex plus</i> *	73 x 89	30	27459
SERVOX® baseplate <i>flex plus</i> *	73 x 89	10	27460
SERVOX® baseplate <i>flex fit</i> *	88 x 108	30	28570
SERVOX® baseplate <i>flex fit</i> *	88 x 108	10	28569

*Due to patent restrictions, these products cannot be used together with the Inhealth HME product BE 1055 in Germany.

7. Disclaimer

Orbisana Healthcare GmbH prohibits any person not authorised by the Orbisana Healthcare GmbH to make any changes to the SERVOX® baseplate or its components. In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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