# SERVOX® Air Pro

## Mesh Nebulizer

DE

EN

FR

CE 0197 MD

REF 38610

38616/Version 2025-04-11

EC REP

Lotus NL B.V.  
Koningin Julianaplein 10,  
1e Verd, 2595AA,  
The Hague, Netherlands

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## Contents

|           |                                                                                            |           |
|-----------|--------------------------------------------------------------------------------------------|-----------|
| <b>1</b>  | <b>Preface .....</b>                                                                       | <b>18</b> |
| 1.1       | Intended purpose.....                                                                      | 19        |
| 1.2       | Intended user.....                                                                         | 19        |
| 1.3       | Indication.....                                                                            | 19        |
| 1.4       | Contraindication.....                                                                      | 19        |
| 1.5       | Side effects .....                                                                         | 19        |
| <b>2</b>  | <b>Safety instructions .....</b>                                                           | <b>19</b> |
| <b>3</b>  | <b>Scope of supply.....</b>                                                                | <b>21</b> |
| <b>4</b>  | <b>Description .....</b>                                                                   | <b>22</b> |
| 4.1       | Charging the battery .....                                                                 | 22        |
| <b>5</b>  | <b>Use .....</b>                                                                           | <b>23</b> |
| 5.1       | Assembly.....                                                                              | 23        |
| 5.2       | Fill inhaler.....                                                                          | 23        |
| 5.3       | Inhalation.....                                                                            | 24        |
| <b>6</b>  | <b>Cleaning and disinfection .....</b>                                                     | <b>25</b> |
| 6.1       | Cleaning.....                                                                              | 25        |
| 6.2       | Disinfection.....                                                                          | 26        |
|           | 6.2.1 Chemical disinfection with 75 % – 95 % ethanol solution.....                         | 26        |
|           | 6.2.2 Thermal disinfection by boiling (not applicable for the masks and the housing) ..... | 26        |
| <b>7</b>  | <b>Control / exchange of consumables.....</b>                                              | <b>26</b> |
| <b>8</b>  | <b>Troubleshooting.....</b>                                                                | <b>27</b> |
| <b>9</b>  | <b>Conditions for operation, storage and transport .....</b>                               | <b>27</b> |
| <b>10</b> | <b>Accessories and spare parts.....</b>                                                    | <b>27</b> |
| <b>11</b> | <b>Useful life/ Disposal .....</b>                                                         | <b>28</b> |
| <b>12</b> | <b>Technical data .....</b>                                                                | <b>28</b> |
| 12.1      | Device specification .....                                                                 | 28        |
| 12.2      | Electromagnetic compatibility (EMC).....                                                   | 29        |
| <b>13</b> | <b>Warranty / Disclaimer .....</b>                                                         | <b>31</b> |

## 1 Preface

This instruction of use is intended for patients, relatives and medical personnel. The SERVVOX® Air Pro Mesh Nebulizer can be used in domestic and clinical settings.

Exclusively the attending physician determines the medication, dosage and frequency of application and the duration of therapy. The SERVVOX® Air Pro Mesh Nebulizer must be used in accordance with the descriptions in this user manual and in accordance with recognized medical regulations.

If you have any questions or problems with the SERVVOX® Air Pro Mesh Nebulizer, please contact your sales representative or Orbisana Healthcare GmbH.

#### **Note**

Read these operating instructions carefully. It contains a lot of important information for a safe and proper application.

For future reference, keep these instructions for use in a safe place near the device.

### **1.1 Intended purpose**

The SERVVOX® Air Pro Mesh Nebulizer nebulizes water-soluble drugs and saline solutions for inhalation for the upper and lower respiratory airways. Pentamidine-containing drugs or medications or liquids that contain oils, essential oils, esters, herbal extracts or dissolved particles or that have corrosive properties must not be nebulized with the SERVVOX® Air Pro Mesh Nebulizer.

The SERVVOX® Air Pro Mesh Nebulizer must not be used for respiration or anesthesia of patients.

EN

#### **Note**

Please only dilute the medication with saline solution.

### **1.2 Intended users and facilities**

The device is a mesh nebuliser designed for the nebulisation of medication for inhalation therapy in professional healthcare facilities and in the home environment. It is suitable for children ( $\geq 3$  years old) and weakened persons and should be used under adult supervision. Before use, instruction in the correct inhalation technique and used by the attending physician or medical personnel is necessary.

#### **Note**

The SERVVOX® Air Pro Mesh Nebulizer is intended for reuse on a single patient and may not be reused by other persons without prior professional reprocessing by the manufacturer or authorized dealers. There is a risk of infection.

### **1.3 Indication**

For respiratory disorders such as asthma and chronic obstructive pulmonary disease (COPD). For colds or other inflammations of the respiratory tract as well as mucus formation.

### **1.4 Contraindication**

The SERVVOX® Air Pro Mesh Nebulizer should not be used by patients with pulmonary edema or by patients with an acute asthma attack or an acute pulmonary infarction.

### **1.5 Side effects**

Various medications can induce a cramping of the respiratory tract (bronchospasm) during inhalation. Ask your attending physician about this.

## **2 Safety instructions**

Safe and proper use of the SERVVOX® Air Pro Mesh Nebulizer is important for your health. If you do not observe these safety instructions, this can have serious health consequences for you and even result in life-threatening injuries.

**Caution:** the word „caution“ indicates a situation, which, if not avoided, may result in minor or moderate injury.

In the following, warnings for safe use are marked with  the symbol.

**Warning:** the word „warning“ alerts you to a situation, which, if not avoided, could result in serious injury or even death.

**Please observe the following safety instructions:**

#### Warning

**Risk of injury:** Children or persons with cognitive disorders may not be able to assess the potential risk of injury when handling the SERVOX® Air Pro Mesh Nebulizer and may injure themselves.

- Do not leave children or persons with cognitive disorders alone with the SERVOX® Air Pro Mesh Nebulizer or the USB charging cable!
- Keep the SERVOX® Air Pro Mesh Nebulizer out of the reach of children or persons with cognitive disorders!

**Risk of electric shock:** If the SERVOX® Air Pro Mesh Nebulizer is damaged or submerged in water or water has penetrated into the interior of the unit, you may receive an electric shock. You can get hurt.

- Do not clean the nebulizer under running water or immerse it in water, especially make sure that the USB port does not get wet!
- Do not use the nebulizer in the shower or when bathing!
- Do not use the nebulizer if it is damaged or liquid has penetrated into it!
- Do not open or modify the device!
- Do not open or replace the battery yourself!
- Only use original spare parts!

**Risk of electric shock:** If the charging adapter is damaged or submerged in water or water has penetrated into the interior of the unit, you may receive an electric shock. You can get hurt.

- Do not clean the adapter under running water or immerse it in water!
- Do not use the adapter if it is damaged or liquid has penetrated into the adapter!
- Do not open or modify the adapter!
- Do not use the adapter if the contacts are dusty or otherwise dirty or damp!

**Fire hazard:** Possible sparks from other electrical appliances or heat sources such as ovens, candles, etc. may cause the aerosol mist to ignite. They can burn themselves or objects can catch fire.

- Do not use the nebulizer near devices that can generate sparks, such as microwave ovens or hair dryers!
- Do not use the nebulizer near open flames, such as candles or heat sources, such as ovens, etc.!
- Do not use the nebulizer in an environment with very high temperatures, e.g. in a sauna!

**Danger of explosion:** If the SERVOX® Air Pro Mesh Nebulizer is exposed to flammable gases, e.g. oxygen or anesthetic mixtures during operation, it may ignite the gases and cause an explosion. You can get hurt.

- Do not use the Nebulizer in areas where flammable gas, oxygen or aerosol products are used!

## Caution

If you wind the USB cable around the nebulizer, the cable may be damaged or broken.

- Do not wrap the USB cable around the nebulizer!

Improper charging may damage or destroy the battery.

- To charge the battery, only use the enclosed charging adapter or chargers with the output value specified by the manufacturer (output 5.0V/1.0A) and compliant with IEC 60601!

Oily substances, esters or aggressive/corrosive substances can stick or destroy the mesh net.

- Do not nebulize oils, herbal extracts, essential oils or aggressive/corrosive substances!

If the nebulizer is exposed to strong vibrations or shocks, the mechanics or the mesh net can be destroyed. The mesh net in the nebulizer unit is sensitive and can be damaged by e. g. direct contact or contact with aggressive cleaning agents or aggressive gases.

- Do not drop the device!
- Protect the device from strong vibrations, e.g. during transport!
- Do not touch the mesh net with your fingers or other objects!
- Do not use aggressive cleaning agents or disinfectants!

EN


## Note

„Please inform us and, if necessary, the competent authority if life-threatening illness or injury occurs in connection with this product or if inpatient treatment results from its use.“

## 3 Scope of supply

## Note

The illustrations are schematic and may therefore differ from the actual view of the SERVOX® Air Pro Mesh Nebulizer.

- 1 Medication cup lid
- 2 Nebulizer unit
- 3 Medication cup
- 4 Nebulizer nozzle with mesh net
- 5 Housing
- 6 Power button
- 7 Switch guard
- 8 Locking button
- 9 Micro-USB-Port
- 10 USB cable
- 11 Mouthpiece, 1 piece
- 12 Mask, 2 piece
- 13 Charging adapter

Also included:  
Case, Instruction of use

Figure 1 Scope of supply

# 4 Description

The SERVOX® Air Pro Mesh Nebulizer nebulizes medicaments into an aerosol, to make them inhalable. A rechargeable battery drives the SERVOX® Air Pro Mesh Nebulizer.

The SERVOX® Air Pro Mesh Nebulizer consists of a hand-held housing and a removable nebulizer unit, consisting of the medication container and the nebulizer nozzle.

The inhalant is filled in the drug container of the nebulizer unit. The nebulizer unit is mounted on the top of the SERVOX® Air Pro Mesh Nebulizer and locked with the locking button on the back of the housing. The inhalate flows down to the nebulizer nozzle. The nebulizer nozzle is equipped with a very finely perforated stainless steel mesh (hole size 3 µm). When the SERVOX® Air Pro is switched on, this mesh network, driven by a piezoelectric ceramic, is set into very fast oscillations. The inhalate is pressed through the fine pores of the mesh net and thus atomized into even small droplets, a so-called aerosol. Aerosols can reach the bronchi and alveoli of the lungs due to their extremely small diameter.


The mouthpiece or mask is mounted on the outside of the nebulizer nozzle of the nebulizer unit. The inhalant is inhaled through the mouthpiece or the mask. A switch guard protects the power button against accidental pressing.

## 4.1 Charging the battery

The SERVOX® Air Pro Mesh Nebulizer is equipped with an internal rechargeable lithium titanate battery. With a fully charged battery, the SERVOX® Air Pro Mesh Nebulizer has a runtime of approx. 60 minutes. At a charging current of 5V/1A the battery is charged after approx. 1 hour.

When the battery is fully charged, the backlight of the power button [6] will glow evenly green during use. When the battery power decreases, the green backlight of the power button starts flashing. At the latest then the battery should be recharged. If the device is stored unused for more than one month, the battery should be charged at least once a month.

To charge the battery, connect the micro-USB plug of the USB cable [10] to the micro-USB port [9] on the back of the housing and the charging adapter [13] to the USB port of the cable and plug it into an electrical outlet.

**Figure 2** Charging the battery

## 5 Use

### Note

- Before first use, clean and disinfect the housing and the individual parts according to the instructions in this instruction of use.
- Charge the battery fully before first use.

### 5.1 Assembly

- Unpack the housing and check for damage.
- Clean and disinfect the housing [5] and all individual parts (see chapter Cleaning and Disinfection).
- Check whether the nebulizer unit [2] is firmly engaged in the housing. If necessary, press the nebulizer unit on the housing until it clicks into place with an audible click.
- Check whether the lid [1] is screwed firmly onto the nebulizer unit. If necessary, tighten with a slight clockwise rotation.
- Push the mouthpiece [11] or mask [12] with light pressure over the nebulizer nozzle [4] of the nebulizer unit


Figure 3 Assembly


### 5.2 Fill inhaler

### Note

- The duration and frequency of use, as well as the dosage and, if necessary, the degree of dilution of the inhalate are determined exclusively by the attending physician. The inhalation solution must be prepared according to the instructions of the attending physician.
- Always hold the nebulizer vertically, when the inhalate is in the medication cup. Do not tilt, rotate or shake the nebulizer too much.
- Discard any Inhalation residues that remain in the medication cup after inhalation. Never reuse them!

- Prepare the inhalation solution according to the attending physician's instructions.
- Check whether the nebulizer unit is firmly engaged in the housing of the Mesh Nebulizer.
- Remove the lid of the medication cup with a slight counterclockwise turn.
- Place the nebulizer vertically and pour a maximum of 10 ml of the inhalant into the medication cup [3]. The medication cup has a scale up to 10 ml. After filling in the inhalant, always hold the mesh nebulizer vertically.
- Tighten the lid of the medication cup with a slight clockwise twist.

EN

**Figure 4** Filling the medication cup

### 5.3 Inhalation

#### Warning


Stop inhalation if sudden discomfort or shortness of breath occurs and consult your attending physician.

#### Note

Always hold the nebulizer upright during inhalation. Do not tilt, rotate or shake the nebulizer excessively. This is the only way to keep the inhaler in contact with the nebulizer unit.

- Sit down with your torso upright.
- **With the mouthpiece:** Take the mouthpiece between the teeth for the duration of inhalation and enclose it tightly with the lips.
- **With the mask:** Keep the mask close to mouth and nose for the duration of inhalation.
- Switch on the Mesh Nebulizer by pressing the power button [6] on the front of the housing and breathe in and out slowly, deeply and evenly. The backlight of the power button will glow evenly green during use.
- If there is no more inhalate in the medication cup, the device switches off automatically. The backlight of the power button winks out when the unit is switched off.

**Figure 5** Inhalation

#### Children and people in need of care:

Children or persons with physical or mental disabilities need the supervision and help of relatives or medical personnel during application.

Use the soft child mask for babies and children. If children scream or cry during inhalation, stop inhalation until quiet and even breathing is possible again.

### Note

The nebulizer switches itself off automatically after 10 minutes of operation. If you want to continue inhalation, make sure that there is still enough inhalate in the medication cup and press the power button again.

## 6 Cleaning and disinfection

### Caution

- Do not use aggressive cleaning agents or disinfectants! Do not use strong oxidants such as perchlorates, or disinfectants that have corrosive effects on metals, macromolecular compounds, or polymers.
- Do not submerge the housing or the charging adapter in water or other liquids.
- Do not disinfect any part of the SERVOX® Air Pro Mesh Nebulizer in the microwave!

EN

Clean the mask, mouthpiece and nebulizer unit after each use and disinfect in the intervals listed below. Normally, the charging adapter does not have to be cleaned or disinfected.

### 6.1 Cleaning

- Switch off the nebulizer.
- Remove the mask or mouthpiece from the nebulizer unit.
- Pour out any residual medication from the medication container and rinse it out with clear tap water.
- Fill the medication cup with some clear tap water, switch on the nebulizer and nebulize all of the water to remove any remaining medication from the mesh.
- Unlock the nebulizer unit by pressing the locking button [8] on the back of the housing and remove it from the housing. Remove the lid of the medication cup.

**Figure 6** Removing the medication cup

- Rinse the nebulizer unit and lid thoroughly with warm tap water (up to 40 °C).
- Soak the mouthpiece or mask for 5 minutes with clean warm tap water (up to 40 °C).
- After cleaning, thoroughly dry the parts with dry sterile medical gauze or a clean, lint-free cloth.
- Wipe the housing with damp sterile medical gauze or a clean, lint-free cloth. Wipe off any inhalate residue on the contact electrodes. Never submerge the housing and in water or put it under running water.
- Dry the housing with dry sterile medical gauze or a clean, lint-free cloth and let it dry completely in the air.

## 6.2 Disinfection

### Note

- If the parts are not cleaned and disinfected correctly and as frequently as indicated, microorganisms may remain in the nebulizing unit, causing risk of infection.

Disinfect the medication container, mouthpiece and lid before using the nebulizer for the first time or after not using the nebulizer for a long period. Clean Mesh Nebulizer as described before disinfection.

### 6.2.1 Chemical disinfection with 75 % – 95 % ethanol solution

- **Nebulizer unit and the mouthpiece or mask:** soak these parts for 5-10 minutes in the ethanol solution or wipe the surfaces of the parts with a lint-free cloth soaked in the ethanol solution.
- Afterwards, rinse the parts with distilled water and then dry them with sterile medical gauze or a clean, lint-free cloth or let them air dry.
- **Housing:** wipe the surface of the housing with a lint-free cloth soaked in the ethanol solution.
- Afterwards dry it with sterile medical gauze or a clean, lint-free cloth or let it air dry.

### 6.2.2 Thermal disinfection by boiling (not applicable for the masks and the housing)

- Due to the material's difference of each part, there are difference on the boiling disinfection. Please refer to the table below:

| Parts                           |                                                                                     | Materials                                             | Boiling               | Time of Boiling | Time interval              |
|---------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------|-----------------------|-----------------|----------------------------|
| Lid                             |    | PC                                                    | Yes                   | ≤10 mins        | Daily during normal use.   |
| Medication cup                  |    | PC and Metal alloy                                    | Yes                   | ≤10 mins        | Monthly during normal use. |
| Mouthpiece under Medication cup |    | PP                                                    | Yes                   | ≤10 mins        | Daily during normal use.   |
| Adult Mask                      |   | PVC                                                   | Do not boil the masks |                 |                            |
| Child Mask                      |                                                                                     |                                                       |                       |                 |                            |
| Housing                         |  | Do not boil the housing, the USB cable or the adapter |                       |                 |                            |

## 7 Control / exchange of consumables

The SERVOX® Air Pro Mesh Nebulizer is generally maintenance-free. Replace the medication cup, the mouthpiece or mask when they show signs of wear after some time due to cleaning and disinfection like, being hard or brittle and blind or if the nebulizer unit is worn out generally.

## 8 Troubleshooting

| Failure                            | Correction                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Does not start when switched on    | <ul style="list-style-type: none"> <li>Check if the battery is charged.</li> <li>Check whether the medication cup is filled sufficiently.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Low fogging performance            | <ul style="list-style-type: none"> <li>Check that the inhalate is in sufficient contact with the mesh of the nebulizer unit: If necessary, tilt the nebulizer slightly forward to allow the inhalate to flow directly to the mesh net.</li> <li>Check if the mesh net is clogged:</li> <li>If the mesh net is blocked, switch off the mesh nebulizer, remove the mouthpiece or mask, dilute 2-3 drops of white household vinegar with 3-6 ml water, pour into the medication cup and completely nebulize (do not inhale) and then clean.</li> <li>Is the inhalate suitable for the Mesh Nebulizer? The inhalate should be water soluble, non-corrosive and free of oils.</li> </ul> |
| Unusual noises during nebulization | <ul style="list-style-type: none"> <li>Check if the filled inhalate is suitable for the Mesh Nebulizer. If necessary, remove inhalant and clean Mesh Nebulizer.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Device stops nebulization          | <ul style="list-style-type: none"> <li>Check to see if there is any inhalate left in the medication cup.</li> <li>Check if the filled inhalate is suitable for the Mesh Nebulizer. If necessary, remove inhalant and clean Mesh Nebulizer.</li> <li>Check if the battery needs charging</li> <li>The mesh nebulizer switches itself off automatically after 10 minutes of operation.</li> </ul>                                                                                                                                                                                                                                                                                     |

If you are unable to correct the failure using the troubleshooting instructions, contact your local sales representative or Orbisana Healthcare GmbH.

EN

## 9 Conditions for operation, storage and transport

Store or transport the SERVOX® Air Pro Mesh Nebulizer protected from direct sunlight, dry and dust-free.

| Characteristic    | Operation                 | Storage and transport     |
|-------------------|---------------------------|---------------------------|
| Temperature       | 5 °C – 40 °C              | - 20 °C – 55 °C           |
| Relative Humidity | ≤ 80 % RF, non condensing | ≤ 80 % RF, non condensing |
| Air pressure      | (86.0 – 106.0) kPa        | (70.0 – 106.0) kPa        |

## 10 Accessories and spare parts

| Articles          | Item number |
|-------------------|-------------|
| mouthpiece        | 38614       |
| Mask for adults   | 38612       |
| Mask for children | 38613       |
| lid               | 38619       |
| nebulizer unit    | 38611       |
| USB cable         | 38615       |
| standard adapter  | 39616       |

Masks, mouthpieces and other accessories or spare parts can be obtained from your sales representative or from Orbisana Healthcare GmbH.

# 11 Useful life/ Disposal

The SERVOX® Air Pro Mesh Nebulizer has an estimated service life of approx. 2 years if used as intended and if all instructions in the instructions for use are observed.

The SERVOX® Air Pro Mesh Nebulizer is subject to Directive 2012/19/EU (WEEE) and must therefore not be disposed of with normal household waste. For proper disposal, the SERVOX® Air Pro Mesh Nebulizer should therefore be cleaned, disinfected, and returned to the manufacturer or an authorized specialist company.

Spare parts such as mouthpieces or masks as well as the packaging can be disposed of with normal household waste.

## 12 Technical data

### 12.1 Device specification

| Specification                                                  |                                                                                                                                                             |
|----------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power supply                                                   | DC 4.8 V (internal lithium titanate battery)<br>DC 5.0 V with mains adapter                                                                                 |
| Power consumption                                              | < 4.0 W                                                                                                                                                     |
| Battery life                                                   | approx. 2 years                                                                                                                                             |
| Charging adapter                                               | <b>Model:</b> LXCP12X.050100BG<br><b>Input:</b> AC 100 – 240 V, 50/60 Hz, 0.5 A Max.<br><b>Output:</b> DC 5 V 1 A USB<br><b>Size:</b> 61.3 x 40.0 x 48.2 mm |
| Atomization rate: according to EN ISO 27427                    | > 0.22 ml/min.                                                                                                                                              |
| Particle Size: according to EN ISO 27427                       | MMAD < 5 µm                                                                                                                                                 |
| Aerosol output according to EN ISO 27427                       | 1.77 mg                                                                                                                                                     |
| Respirable fraction (%; 0.5 – 5 µm): according to EN ISO 27427 | 72 % ± 5 %                                                                                                                                                  |
| Operating frequency                                            | 110 kHz ± 10 kHz                                                                                                                                            |
| Operating noise                                                | 40 dB ± 5 dB                                                                                                                                                |
| Permissible temperature of the inhalate                        | ≤ 45 °C                                                                                                                                                     |
| Capacity of the medicine cup                                   | 10 ml                                                                                                                                                       |
| Size                                                           | 37.6 mm (L) × 37.6 mm (W) × 97.7 mm (H)                                                                                                                     |
| Weight                                                         | 100 g (incl. Akku)                                                                                                                                          |
| Security level                                                 | Internal power supply Type BF                                                                                                                               |
| Working Environment                                            | <b>Temperature:</b> 5 °C – 40 °C<br><b>Relative Humidity:</b> ≤ 80 % R.H. Non-condensing state<br><b>Atmospheric pressure:</b> (86.0 ~ 106.0) kPa           |
| Storage/delivery Environment                                   | <b>Temperature:</b> - 20 °C – 55 °C<br><b>Relative Humidity:</b> ≤ 80 % R.H. Non-condensing state<br><b>Atmospheric pressure:</b> (70.0 ~ 106.0) kPa        |

## 12.2 Electromagnetic compatibility (EMC)

### Guidance and Manufacturer's declaration – electromagnetic EMISSIONS

This device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such environment.

#### Electromagnetic emission IEC60601-1-2

| Emissions text                                        | Compliance | Electromagnetic environment-guidance                                                                                                                                                                                     |
|-------------------------------------------------------|------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions CISPR11                                  | Group 1    | This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                         |
| RF emissions CISPR11                                  | Class B    | This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                      | Class A    |                                                                                                                                                                                                                          |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Complies   |                                                                                                                                                                                                                          |

EN

### Guidance and Manufacturer's declaration – electromagnetic IMMUNITY

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such an environment.

| Immunity test                                                                           | IEC 60601 test level                                                                                                                                                | Compliance level                                                                                                                                                      | Electromagnetic environment-guidance.                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD) IEC 61000-4-2                                             | ± 6 kV contact<br>± 8 kV air                                                                                                                                        | ± 6 kV contact<br>± 8 kV air                                                                                                                                          | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                  |
| Electrical fast transient/burst IEC 61000-4-4                                           | ± 2kV for power supply lines<br>± 1 kV for input/output lines                                                                                                       | ± 2 kV for power supply lines*1)                                                                                                                                      | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                             |
| Surge IEC 61000-4-5                                                                     | ± 1 kV line to line<br>± 2 kV line to earth                                                                                                                         | ± 1 kV line to line<br>± 2 kV line to earth                                                                                                                           | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                             |
| Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11 | < 5 % UT (> 95 % dip in UT) for 0.5 wde<br>40 % UT (60 % dip in UT) for 5 cycles<br>70 % UT (30 % dip in UT) for 25 cycles<br>< 5 % UT (> 95 % dip in UT) for 5 sec | < 5 % UT (> 95 % dip in UT) for 0.5 cycle<br>40 % UT (60 % dip in UT) for 5 cycles<br>70 % UT (30 % dip in UT) for 25 cycles<br>< 5 % UT (> 95 % dip in UT) for 5 sec | Mains power quality should be that of a typical commercial and/or hospital environment. If the user of this device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or battery. |

## Guidance and Manufacturer's declaration – electromagnetic IMMUNITY

The device is intended for use in the electromagnetic environment specified below.  
The customer or the user of the device should assure that it is used in such an environment.

| Immunity test                                              | IEC 60601 test level        | Compliance level | Electromagnetic environment-guidance.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------------------------------------------|-----------------------------|------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8 | 3 A/m                       | 3 A/m            | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. Portable and mobile RF communications equipment should be used no closer to any part of this device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.                                                                                                                                                                            |
| Conducted RF<br>IEC 61000-4-6                              | 3 V rms<br>15 KHz to 80 MHz | 3 V/m            | <b>Recommended separation distance:</b><br>$d = 1.2 \sqrt{P}$<br>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Radiated RF<br>IEC 61000-4-3                               | 3 V/m<br>80 MHz to 2.5 GHz  | 3 V/m            | Recommended 80 MHz to 800 MHz separation distance<br>$d = 800 \text{ MHz to } 2, 5 \text{ GHz where } P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in meters (m).<br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,*2) should be less than the compliance level in each frequency range.*3) Interference may occur in the vicinity of ((( ))) equipment marked with the following symbol: ((( ))) |

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating this device.

- Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. Recommended separation distances between portable and mobile RF communications equipment and the device.  
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

| maximum output power rate of transmitters (w) | separation distance according to frequency of transmitter |                                |                                |
|-----------------------------------------------|-----------------------------------------------------------|--------------------------------|--------------------------------|
|                                               | 150 kHz to 80 MHz<br>$d = 1.2$                            | 80 MHz to 800 MHz<br>$d = 1.2$ | 800 MHz to 2.5 GHz<br>$d = 23$ |
| 0.01                                          | 0.12                                                      | 0.12                           | 0.23                           |
| 0.1                                           | 0.38                                                      | 0.38                           | 0.73                           |
| 1                                             | 1.2                                                       | 1.2                            | 2.3                            |
| 10                                            | 3.8                                                       | 3.8                            | 7.3                            |
| 100                                           | 12                                                        | 12                             | 23                             |

EN

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

#### Caution

- In order to regulate the requirements for EMC with the aim to prevent unsafe product situations, the EN60601-1-2 standard has been implemented. The nebulizer Air Pro conforms to the EN60601-1-2:2012 standard for both immunity and emissions.
- Portable and mobile RF communication equipment may affect the performance of the nebulizer; avoid strong electromagnetic interference when using, such as close to mobile phones, microwave ovens and so on.

## 13 Warranty / Disclaimer

The manufacturer's warranty covers manufacturing and material defects and is valid for one year from the date of purchase for the main unit, six months from date of purchase for the medicine cup and three months from the date of purchase for the remaining parts and accessories. Defective parts, with the exception of wear and spare parts are either replaced or repaired by the manufacturer. It is for the manufacturer to decide. If you have any questions or problems with the product, please contact your provider or contact person at Orbisana Healthcare GmbH.

The warranty claim expires if the care and maintenance instructions described in these operating instructions are not observed. The manufacturer forbids any modifications to this device, otherwise the warranty claim will also expire.

In the event of non-observance of the instructions for use or failure to observe the safety instructions or manipulation of the product, the manufacturer assumes no liability in the event of damage.

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