SERVOX® | SERVONA

Absaugkatheter

suction catheter

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Catéter de succión

usisni kateter

Sací katéter

DE

EN

FR

IT

ES

HR

SK

39882/Version 2024-11-19

CE 0197 MD


Orbisana Healthcare GmbH

Biberweg 24-26
53842 Troisdorf | Germany
Tel. +49 (0)2241 9322-0
Fax +49 (0)2241 9322-271
info@orbisana.de
www.orbisana.de

DE: SERVONA Absaugkatheter	ES: SERVONA Catéter de aspiración
EN: SERVONA Suction Catheter	HR: SERVONA usisni kateter
FR: SERVONA Cathéter d'aspiration	SK: Sací katéter SERVONA
IT: SERVONA Catetere aspirazione	

15 cm		Shore	Chariere (CH)	REF
	15 cm / 11 cm	75	CH 6	39741
	15 cm / 11 cm	75	CH 8	25380
	15 cm / 11 cm	75	CH 10	25381
	15 cm / 11 cm	75	CH 12	39742
	15 cm / 11 cm	75	CH 14	39743
	15 cm / 11 cm	75	CH 16	39744

50 cm		Shore	Chariere (CH)	REF
	50 cm / 46 cm	75	CH 6	22375
	50 cm / 46 cm	75	CH 8	22376
	50 cm / 46 cm	75	CH 10	22377
	50 cm / 46 cm	75	CH 12	22378
	50 cm / 46 cm	75	CH 14	22379
	50 cm / 46 cm	75	CH 16	22380
	50 cm / 46 cm	75	CH 18	22381

DE: SERVVOX® Absaugkatheter	ES: SERVVOX® Catéter de aspiración
EN: SERVVOX® Suction Catheter	HR: SERVVOX® usisni kateter
FR: SERVVOX® Cathéter d'aspiration	SK: Sací katéter SERVVOX®
IT: SERVVOX® Catetere aspirazione	

30 cm soft		Shore	Chariere (CH)	REF
	30 cm / 26 cm	65	CH 6	39745
	30 cm / 26 cm	65	CH 8	39746
	30 cm / 26 cm	65	CH 10	39747
	30 cm / 26 cm	65	CH 12	39748
	30 cm / 26 cm	65	CH 14	39749
	30 cm / 26 cm	65	CH 16	39750
	30 cm / 26 cm	65	CH 18	39751

DE: SERVOX® Absaugkatheter FT	ES: SERVOX® Catéter de aspiración FT
EN: SERVOX® Suction Catheter FT	HR: SERVOX® usisni kateter FT
FR: SERVOX® Cathéter d'aspiration FT	SK: Sací katéter SERVOX® FT
IT: SERVOX® Catetere aspirazione FT	

50 cm		Shore	Chariere (CH)	REF
	50 cm / 43 cm	75	CH 6	39752
	50 cm / 43 cm	75	CH 8	39753
	50 cm / 43 cm	75	CH 10	39754
	50 cm / 43 cm	75	CH 12	39755
	50 cm / 43 cm	75	CH 14	39756
	50 cm / 43 cm	75	CH 16	39757
	50 cm / 43 cm	75	CH 18	39758

DE

DE: SERVOX® Absaugkatheter mit atraumatischer Spitze FT	ES: SERVOX® Catéter de aspiración con punta atraumática FT
EN: SERVOX® Suction Catheter with atraumatic tip FT	HR: SERVOX® usisni kateter s atraumaticnim vrhom FT
FR: SERVOX® Cathéter d'aspiration à pointe atraumatique FT	SK: Sací katéter SERVOX® s atraumatickým hrotom FT
IT: SERVOX® Catetere aspirazione con punta atraumatica FT	

50 cm		Shore	Chariere (CH)	REF
	50 cm / 43 cm	75	CH 8	39760
	50 cm / 43 cm	75	CH 10	39761
	50 cm / 43 cm	75	CH 12	39762
	50 cm / 43 cm	75	CH 14	39763
	50 cm / 43 cm	75	CH 16	39764

DE: SERVOX® Absaugkatheter atraumatisch mit Trichteranschluss	ES: Catéter de succión SERVOX® atraumático con conector de embudo
EN: SERVOX® suction catheter atraumatic with funnel connection	HR: SERVOX® usisni kateter atraumatski s priključkom za lijevak
FR: SERVOX® Cathéter d'aspiration atraumatique avec raccord entonnoir	SK: Sací katéter SERVOX® atraumatický s lieviovým pripojením
IT: SERVOX® Catetere di aspirazione atraumatico con raccordo a imbuto	

50 cm		Shore	Chariere (CH)	REF
	50 cm / 46 cm	75	CH 8	26264
	50 cm / 46 cm	75	CH 10	26265
	50 cm / 46 cm	75	CH 12	26266
	50 cm / 46 cm	75	CH 14	26267
	50 cm / 46 cm	75	CH 16	26268

39882/Version 2024-11-19

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SERVOX® | SERVONA Absaugkatheter

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EN

1 Foreword

This instruction manual is intended for patients, relatives and medical staff. The patient can use the SERVOX® | SERVONA suction catheter independently (without the support of relatives or carers), but must be given instructions on its use by their doctor or by medical staff beforehand. SERVOX® | SERVONA suction catheters can be used at home and in a clinical setting. They are suitable for all age groups according to the size selected by the patient's doctor.

SERVOX® | SERVONA suction catheters should be used in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine.

Note

Read this instruction manual carefully. It contains important information on safe and proper use of the product.

39882/Version 2024-11-19

SERVOX® | SERVONA suction catheters are assigned to batches for traceability purposes. The batch number of each suction catheter is indicated on the packaging label next to "LOT". This batch number allows each suction catheter to be identified.

If you have any questions or problems with the suction catheters, please contact your field representative or Orbisana Healthcare GmbH, referencing this batch number.

1.1 Intended purpose

The SERVOX® | SERVONA suction catheters are intended for suctioning secretions from the respiratory tract and trachea. They are suitable for nasal, oral and endotracheal suctioning.

Note

Suction catheters are sterile disposable products and must not be re-sterilised or reused.

1.2 Indication

For use on ventilated or spontaneously breathing tracheotomised, laryngectomised or intubated patients or patients who are not tracheotomised, laryngectomised or intubated but require help with expectoration.

1.3 Contraindication

Do not use on patients with skull fractures or patients with known allergies to the material used.

Patients with an acute tendency to bleed due to e.g. anticoagulation treatment may only use the SERVOX® | SERVONA suction catheter after consulting their doctor.

1.4 Side effects

There are no known side effects provided the patient uses the product properly.

2 Safety instructions

Safe and proper use of your SERVOX® | SERVONA suction catheter is important for your health. Ignoring the safety instructions can have significant effects on your health or even cause life-threatening injuries.

In the following, warnings on the safe use of the product are indicated with the  symbol.

Danger: The word "Danger" indicates a likely situation that, if not avoided, can result in grave injury or even death.

Warning: The word "Warning" indicates a situation that, if not avoided, can result in grave injury or even death.

Please observe the following safety instructions:

Danger

Risk of injury/ suffocation: The tube adhering to the tracheal wall, careless suctioning, suctioning too deep or using damaged suction catheters can irritate or injure the mucosa. This can cause the mucosa to swell and the trachea to seal, which can lead to suffocation. The vagus nerve can be irritated. This can cause cardiac arrest.

- Never use a damaged suction catheter!
- Always check the suction depth during suctioning!
- Always check the suction strength during suctioning!
- Keep the suction catheter moving during suctioning!
- During suctioning, try not to touch the trachea!

Risk of infection: The use of dirty suction catheters or unhygienic work practices can cause infections.

- Do not use the suction catheter after the expiry date!
- Do not use the suction catheter if the packaging is damaged!
- Do not use a used suction catheter!
- Always ensure that you follow hygienic work practices!
- Wear disposable gloves during suctioning!

Danger of suffocation: Selecting the incorrect diameter of the suction catheter can interrupt breathing or cause negative pressure in the lungs.

- Always check the size of the suction catheter!

EN

Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.



3 Description

SERVOX® | SERVONA suction catheters are free of latex and phthalate and are available in various lengths, sizes, hardness grades and models. All SERVOX® | SERVONA suction catheters have a rounded tip.

SERVOX® | SERVONA suction catheters are used together with a suction unit and are comprised of a flexible PVC tube with a connection piece to connect it to the suction unit. To check the insertion depth, a scale in 1-cm increments is marked on the tube.

The lateral eyes on the sides of the tube above the rounded distal tip of the suction catheter prevent it from adhering to the mucosa. The atraumatic tip also has a bulbous head on the distal end of the catheter in order to minimise contact with the mucosa. The various tip shapes differ depending on how they are inserted into the trachea.

Connections:

	
Funnel: The funnel connection must be connected to the suction tube of a suction unit with an additional fingertip control.	Fingertip: The fingertip connection can be used to connect the suction catheter to the suction tube on the suction unit.

Tips:

	
Atraumatic tip: Bulbous head, rounded edge, 4 small offset eyes	Standard tip: Rounded edge, 2 large offset eyes or 2 small eyes opposite each other
➔ Insert with suction	➔ Insert without suction

Size indicated by the coloured marking on the connection pieces:

Size in Charrière (CH)	CH 6	CH 8	CH 10	CH 12	CH 14	CH 16	CH 18
Outer diameter in mm	2.0	2.67	3.33	4.0	4.67	5.33	6.0
Colour	Light green	Light blue	Black	White	Green	Orange	Red

4 Use

Note

- Always insert the suction catheter with the standard tip without suction.
- Always insert the suction catheter with the atraumatic tip with suction.

- Disinfect hands and put on gloves
- Check that the size, length and model of the suction catheter is correct
- Check the sterile packaging for the expiry date and damage. Do not use if the suction catheter has expired or the packaging is damaged
- Open the sterile packaging, take out the connection piece and slide it on the adapter or the fingertip control of the suction tube
- Switch on the suction unit (read the suction unit instruction manual)
- Then remove the rest of the suction catheter from the packaging

- Insert the suction catheter with the **standard tip** carefully **without suction!**
- Insert the suction catheter with the **atraumatic tip** carefully **with suction!**
- Then slowly pull the catheter out under suction with a slight twisting motion
- Do not suction for more than 15 s
- Dispose of the suction catheter after use

5 Storage, transport, disposal

Store or transport the SERVUX® | SERVONA suction catheters in a clean, dry environment that is free of dust and protected from light at temperatures between - 5 °C and + 35 °C.

Protect from direct sunlight and heat.

The SERVUX® | SERVONA suction catheters must be disposed of in line with national regulations for waste products. Suction catheters can generally be disposed of in household waste. The packaging can be recycled.

EN

6 Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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