# SERVOX<sup>®</sup> Tracheo Nose Neo F | Neo P

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## Contents

<b>1 Foreword</b> .....	<b>8</b>
1.1 Intended purpose .....	9
1.2 Indications .....	9
1.3 Contraindications .....	9
<b>2 Safe use</b> .....	<b>9</b>
<b>3 Description</b> .....	<b>11</b>
<b>4 Application</b> .....	<b>11</b>
<b>5 Technical data</b> .....	<b>12</b>
<b>6 Storage and transport conditions</b> .....	<b>13</b>
<b>7 Order information</b> .....	<b>13</b>
<b>8 Disclaimer</b> .....	<b>13</b>

## 1 Foreword

This instruction for use are intended for patients, relatives over 16 years of age and trained medical staff. The SERVVOX® Tracheo Nose Neo and the SERVVOX® Tracheo Nose Neo P can be used at home and in a clinical setting. Prior instruction on use by trained medical staff is required. Application should be in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine.

### Note

Read this instruction manual carefully. It contains important information on the safe and proper use of the product. Keep this instruction manual for future reference.

SERVOX® Tracheo Nose Neo F and Neo P are assigned a batch number for the purpose of unambiguous traceability. The batch number is indicated on the packaging next to "LOT".

If you have any questions or problems with the SERVOX® Tracheo Nose Neo F and Neo P, please contact your field representative or Orbisana Healthcare GmbH, quoting the product's batch number.

### 1.1 Intended purpose

The SERVOX® Tracheo Nose Neo F and Neo P are intended for the one-time use in spontaneously breathing intubated or tracheotomised neonates. They warm and moisten the inhaled air and prevent the ingress of dust and foreign bodies. The product is intended for single use only in order to protect users and patients from cross infections.

#### Note

The SERVOX® Tracheo Nose Neo F and the SERVOX® Tracheo Nose Neo P are non-sterile single-use products that must not be reused or reprocessed because of the risk of infection.

EN

### 1.2 Indication

Use in spontaneously breathing intubated or tracheotomised neonates with a tidal volume of 15 – 50 ml.

### 1.3 Contraindication

- Use outside of the recommended tidal volume range of 15-50 ml
- Severe exsiccosis (dehydration)
- Excessive secretion inside the airways
- Concomitant use of a speech valve
- Allergic reaction to one of the components

## 2 Safe use

Safe and proper handling is important for the patient's health! Ignoring the safety in-structions can have significant effects on your health or even cause life-threatening situations.

**Warning:** The word "Warning" indicates a situation that, if not avoided, can result in grave injury or even death.

## Please observe the following safety instructions:

### Warning

**Risk of infection:** Using a contaminated SERVOX® Tracheo Nose Neo can result in an infection.

- Exchange the SERVOX® Tracheo Nose Neo Module in case of heavy mucus secretion or soiling of the filter medium, but latest after 24 hours!
- Do not use the SERVOX® Tracheo Nose Neo if it or the packaging is damaged!
- Do not use a SERVOX® Tracheo Nose Neo that has been dropped but replace it with a new SERVOX® Tracheo Nose Neo!
- Do not use a SERVOX® Tracheo Nose Neo that has been dropped accidentally but replace it with a new SERVOX® Tracheo Nose Neo!
- Never use water for rinsing the SERVOX® Tracheo Nose Neo!
- Dispose of a used SERVOX® Tracheo Nose Neo immediately after removal!

**Risk of suffocation/risk of respiratory distress:** Accidentally or permanently blocking the tracheostoma can result in breathlessness or suffocation of the patient.

- Never insert the SERVOX® Tracheo Nose Neo directly into a tracheostoma!
- Never modify the SERVOX® Tracheo Nose Neo in any way!
- Immediately remove the SERVOX® Tracheo Nose Neo if the filter is clogged or in case of increased airway resistance e.g. due to increased accumulation of secretions!
- Never cover or block the filter opening of the SERVOX® Tracheo Nose Neo with clothing or similar!
- Do not use any atomisers or humidifiers and do not add moisture to the SERVOX® Tracheo Nose Neo with water, as this can increase breathing resistance and block respiration!
- Do not use the SERVOX® Tracheo Nose Neo in water, the shower or the bathtub!

**Risk of allergy:** Remove the SERVOX® Nose Neo immediately at the first signs of an allergic reaction, such as a skin rash, redness, breathlessness or rapid heartbeat!

**Increased CO<sub>2</sub> concentration in the respiratory air:** A tidal air volume range above 15 – 50 ml can result in insufficient moistening of the respiratory air, a tidal air volume range below 15 – 50 ml can result in an increase of CO<sub>2</sub> concentration in the respiratory air.

- Observe the approved tidal air volume range!

## Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.

## 3 Description

The SERVOX® Tracheo Nose Neo is a heat and moisture exchanger (HME) that can be connected to a 15 mm standard connector of a tracheal cannula or an endotracheal tube. Moisture and heat of the expired air are preserved in a filter medium upon exhale and then returned to the breathing air on inhaled. The humidifying effect sets in after just a few breaths. The filter also protects against the ingress of dust, dirt, insects or other foreign bodies with the inhaled air. Using the SERVOX® Tracheo Nose Neo increases the airway resistance.

To enrich the respiratory air with additional oxygen (O<sub>2</sub>), the SERVOX® Tracheo Nose Neo O<sub>2</sub> Adapter can be fitted onto the SERVOX® Tracheo Nose Neo.

The housing of the SERVOX® Tracheo Nose Neo is made from latex- and DEHP-free ABS plastic without metallic or electrically conductive materials. There are no restrictions for use during magnetic resonance tomography (MRT). The filter core in the front of the housing consists of PU foam for the SERVOX® Tracheo Nose Neo F and paper for the SERVOX® Tracheo Nose Neo P.

In a sudden rise of pressure, e.g. as the result of coughing, the filter is blown out to protect the airways. In that case, the SERVOX® Tracheo Nose Neo needs to be replaced.

SERVOX® Tracheo Nose Neo F

SERVOX® Tracheo Nose Neo P

## 4 Application

Wash the hands thoroughly before use.

**Check prior to use:** Ensure that the packaging and the SERVOX® Tracheo Nose Neo are without damage. Do not use if the packaging and the SERVOX® Tracheo Nose Neo are damaged.

EN

**Fitting:** Use light pressure to press the Die SERVOX® Tracheo Nose Neo with the 15 mm connector onto the standard 15 mm connector of the tracheal cannula or the endotracheal tube until it is firmly fitted.

In a sudden rise of pressure, e.g. as the result of coughing, the filter is blown out to protect the airways. In that case, the SERVOX® Tracheo Nose Neo needs to be replaced.

**Removal:** Remove the SERVOX® Tracheo Nose Neo latest after 24 hours and immediately if the filter is clogged or in case of increased airway resistance. To that end, hold onto the 15 mm standard connector of the tracheal cannula or the endotracheal tube and, pulling lightly, disconnect the SERVOX® Tracheo Nose Neo carefully from the 15 mm standard connector and discard immediately as residual waste.

#### Note

If you have any questions or experience complications, please contact your doctor.

## 5 Technical data

	SERVOX® Tracheo Nose Neo F	SERVOX® Tracheo Nose Neo P
Tidal volume	15 – 50 ml	15 – 50 ml
Connections	15 mm standard connector	15 mm standard connector
Moisture loss after 1 hour	7.6 mg/l (for VT = 250 ml) 8.7 mg/l (for VT = 500 ml)	12.5 mg/l (for VT = 250 ml) 13.4 mg/l (for VT = 500 ml)
Pressure drop after 0 hours	2.96 hPa (for 30 l/min) 11.35 hPa (for 60 l/min) 24.32 hPa (for 90 l/min)	2.4 hPa (for 30 l/min) 6.9 hPa (for 60 l/min) 13.2 hPa (for 90 l/min)
Pressure drop after 24 hours	3.46 hPa (for 30 l/min) 13.95 hPa (for 60 l/min) 25.57 hPa (for 90 l/min)	2.9 hPa (for 30 l/min) 8.0 hPa (for 60 l/min) 15.8 hPa (for 90 l/min)
Interior volume:	3 ml	3 ml
Weight	3.3 g	3.3 g
Maximum period of use	24 hours	24 hours

From the date of manufacture, the SERVOX® Tracheo Nose Neo has a shelf life of 5 years if the package is unopened and stored correctly.

## 6 Storage and transport conditions

Protect from moisture, dust and direct sunlight. and heat.

	Storage and transport conditions
Temperature	+5 °C – +25 °C
relative humidity	≤ 65 % RH

## 7 Order information

Designation	PU	Article no.
SERVOX® Tracheo Nose Neo P	30	30027
SERVOX® Tracheo Nose Neo F	30	28391
SERVOX® Tracheo Nose Neo O <sub>2</sub> Adapter	4	29660

EN

## 8 Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product. w

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