# SERVOX<sup>®</sup> Tracheo Nose P

(Typ: BY-4)

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# INSTRUCTIONS FOR USE

## SERVOX® Tracheo Nose P (Typ: BY-4)

Artificial Nose  
Heat and Moisture Exchanger/ HME

Symbol explanation .....

8

EN

## Foreword

These instructions for use are intended for patients, relatives and healthcare professionals. The SERVOX® Tracheo Nose P (Typ: BY-4) can be used in domestic and clinical settings. It should be used as described in these instructions for use and in accordance with recognised rules of medicine.

### Note

Read these instructions for use carefully. Keep it in a safe place for future reference.

The SERVOX® Tracheo Nose P (Typ: BY-4) is assigned a batch code for unambiguous traceability. The batch code is indicated on the packaging behind the marking "LOT". If you have any questions or problems with the SERVOX® Tracheo Nose P (Typ: BY-4), please contact Orbisana Healthcare GmbH, stating the batch number.

### Intended use

The SERVOX® Tracheo Nose P (Typ: BY-4) warms, humidifies and filters the inhaled air by recovering moisture and heat from the exhaled air by means of a filter, thus reducing irritation and illnesses of the respiratory tract.

### Note

The SERVOX® Tracheo Nose P (Typ: BY-4) is a single use device and must not be reused or reprocessed, there is a risk of infection

### Indication

Use in spontaneously breathing intubated or tracheotomised patients with a respiratory volume of 250 – 1.000 ml.

### Contraindications

Severe exsiccosis (dehydration) with excessive secretion of the respiratory tract

41713/Version 2024-11-14

# Safety instructions

Non-observation of these safety instructions can result in serious health effects or even life-threatening situations.

## Warning

**Risk of infection:** If you use a dirty SERVOX® Tracheo Nose P (Typ: BY-4), you can get an infection.

- If the filter medium becomes heavily clogged or soiled, replace it immediately or after 24 hours at the latest!
- Do not use if damaged or if the packaging is damaged and after the expiry date!
- Do not use a HME that have fallen down, but replace them with a new HME!
- Never put accidentally ejected filters (e. g. due to heavy coughing) back into the housing, but use a new HME!
- Do not rinse the HME with water!
- Dispose of used HME immediately after use!

**Risk of suffocation/danger of respiratory distress:** If your tracheostoma is blocked unintentionally or permanently, you may experience respiratory distress or suffocate.

- Never insert the HME directly into a tracheostoma!
- Never modify the HME!
- Never stick the HME to a tracheostomy tube with adhesive tape, as this makes removal more difficult in case of respiratory distress or danger of suffocation.
- Remove the HME immediately if the filter is clogged or if airway resistance increases, e. g. due to increased secretion accumulation!
- Never cover or block the filter opening of the HME with clothing or the like!
- Do not use nebulisers or humidifiers or moisten the HME with water, as this may increase the airway resistance and block breathing!
- Never wear the HME in water, under the shower or in the bathtub!
- Change the HME regularly when using it in humid environments, e. g. in the rain!

**Increased CO<sub>2</sub> concentration in the respiratory air:** RA tidal volume range above or below 250 – 1000 ml can lead to an increase in CO<sub>2</sub> concentration or insufficient humidification.

- Always observe the prescribed breathing volume!

## Note

Please inform us and, if necessary, the competent authority if life-threatening illness or injury occurs in connection with this product or if in-patient treatment results from its use.

## Description

The SERVOX® Tracheo Nose P (Typ: BY-4) is a heat and moisture exchanger (HME) that can be attached to a 15 mm standard connector of a tracheal cannula or endotracheal tube. When the patient exhales, moisture and heat from the breathing air are retained in the filter medium and released back into the breathing air the next time the patient inhales. After a few breaths, the humidification effect sets in. The filter also protects against the penetration of dust, dirt, insects and other foreign bodies with the breathable air. During use, the airway resistance increases, similar to the physiological resistance of the upper airways, thus supporting the natural lung function.

The housing is made of latex and DEHP-free polypropylene (PP) with no metallic or electrically conductive materials, there is no restriction on use during magnetic resonance imaging (MRI). The two filter cores contained in the front of the housing are made of corrugated paper. In the event of a sudden increase in pressure, e. g. due to coughing, the filters are blown out to protect against overpressure in the airways.

# Use

- Wash hands thoroughly.
- Check that the packaging and HME are undamaged.
- **Attach:** Attach the 15 mm connector to the 15 mm standard connector of the tracheal cannula or endotracheal tube firmly enough with light pressure.
- **Remove:** Hold the 15 mm standard connector of the tracheal cannula or endotracheal tube and carefully pull the HME off the 15 mm standard connector and dispose of it immediately in the residual waste.
- **Duration of use:** Replace after 24 hours at the latest; replace immediately in the event of increased mucus build-up or contamination or if breathing resistance increases

## Note

Contact your attending physician if you have any questions or complications

# Technical data

The SERVOX® Tracheo Nose P (Typ: BY-4) is a sterile, disposable medical device.

Internal volume of the HME	12 ml	Recommended tidal volume	250 – 1.000 ml
Moisture loss (1 hour)	15.6 mg/L (VT 250 ml)	Moisture loss (24 hours)	12.8 mg/L (VT 250 ml)
	15.3 mg/L (VT 500 ml)		15.0 mg/L (VT 500 ml)
	17.2 mg/L (VT 750 ml)		16.4 mg/L (VT 750 ml)
	17.8 mg/L (VT 1.000 ml)		15.8 mg/L (VT 1.000 ml)
Pressure drop (Before use)	0.1 hPa (at 30 L/min.)	Pressure drop (After 24 hours)	0.1 hPa (at 30 L/min.)
	0.2 hPa (at 60 L/min.)		0.2 hPa (at 60 L/min.)
	0.1 hPa (at 90 L/min.)		0.4 hPa (at 90 L/min.)
Connections	15 mm standard connector	Condensation surface	86 cm <sup>2</sup> (+ / - 8 cm <sup>2</sup> )

EN

# Storage and transport conditions

Store or transport in a clean, dry environment at a relative humidity less than 80%. Protect from moisture, corrosive gases, and direct sunlight. Can be used unopened for 3 years from the date of manufacture.

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