# SERVOX® Doppel-Trachealkomresse AL

## standard, gelocht

**Orbisana**  
GESUNDHEITSWELT

90 x 100 mm

EN	SERVOX® Double Tracheal Compress AL, standard, with punched hole
FR	SERVOX® Double compresses trachéales AL, standard, perforée
IT	SERVOX® Doppia compressa tracheale AL, standard, perforata
ES	SERVOX® Compresas traqueales dobles AL, estándar, perforadas
BG	SERVOX® Двоен трахеален компрес AL, стандартен, с отвор
HU	SERVOX® standard dupla légszű-sebpráma AL, lyukas
EL	Επίθερα τραχειοστομίας SERVOX® AL, με σχισμή
HR	Dvostruka trahealna kompresa SERVOX®, AL, standardna, perforirana
NL	SERVOX® Dubbele trachealkompres AL standaard, geperforeerd
PT	Compressa traqueal dupla SERVOX® AL standard, perfurada
PL	Podwójny kompres tracheostomijny SERVOX® AL standard, z otworami
SL	Dvojna standardna trahealna kompresa z luknjo SERVOX® AL

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**EN**

**Intended purpose** SERVOX® Double Tracheal Compress AL standard, with punched hole provide cushioning for the neck flange and absorb moisture that secretes from the tracheostoma to prevent skin irritation and infections.

**Indication** Tracheotomised or laryngectomised tracheal cannula users with heavy secretions.

**Description** The aluminium coating helps prevent the compress from sticking to the skin. The punched hole ensures a more secure fit and prevents slipping.

**Application** Slide the tracheal cannula through the punched hole of the compress (with the aluminium coating to the tracheostoma) until the neck flange is well cushioned. Then insert the tracheal cannula as usual.

To remove, pull out the tracheal cannula as usual, pull the compress off the cannula tube and dispose of it. – **Changing frequency:** with every cannula pull more frequently as needed, at least every 24 hours. – If the skin is irritated, use a compress with no aluminium coating and consult a doctor.

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